NEURAMIS®

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COSMETIC

A Randomized Clinical Trial to Evaluate the Efficacy and Safety of Lidocaine-Containing Monophasic Hyaluronic Acid Filler for Nasolabial Folds

Hong Jin Joo, M.D.,
Young Jun Woo, M.D.,
Jung Eun Kim, M.D., Ph.D.,
Beom Joon Kim, M.D., Ph.D.,
Hoon Kang, M.D., Ph.D.

Background: Hyaluronic acid dermal fillers are most frequently used for unwanted wrinkles. Recently, lidocaine has been incorporated into hyaluronic acid fillers to reduce injection discomfort.

Methods: A randomized, multicenter, double-blind, intraindividual trial was designed to compare a new lidocaine-containing monophasic hyaluronic acid filler (Neuramis Deep Lidocaine) with a lidocaine-containing biphasic hyaluronic acid filler (Restylane Perlane-L) in moderate to severe nasolabial folds. Fifty-eight patients with moderate to severe nasolabial folds were randomized to an injection of Neuramis or Perlane-L in the left or right side of the face. Clinical efficacy and safety were assessed by blinded investigators, independent expert panels, and patients based on the Wrinkle Severity Rating Scale and the Global Aesthetic Improvement Scale at weeks 8, 16, and 24 after the injection.

Results: Wrinkle Severity Rating Scale improvement from baseline with Neuramis (1.64 ± 0.74) was significantly greater than with Perlane-L (1.45 ± 0.54) at week 24 ($p < 0.05$). The mean Global Aesthetic Improvement Scale score at week 24 was 2.36 ± 0.55 for Neuramis and 2.00 ± 0.50 for Perlane-L ($p < 0.05$). However, the difference in pain reduction between Neuramis and Perlane-L-treated sides was not statistically significant.

Conclusions: The efficacy and safety of Neuramis are comparable to those of Perlane-L in Wrinkle Severity Rating Scale and Global Aesthetic Improvement Scale improvement for the management of nasolabial folds. Furthermore, the difference in pain reduction between the two fillers was not clinically significant. (*Plast. Reconstr. Surg.* 137: 799, 2016.)

CLINICAL QUESTION/LEVEL OF EVIDENCE: Therapeutic, I.

The aging process of the skin involves the loss of hyaluronic acid, subcutaneous tissue, and collagen and elastic fibers; the definitive sign of aging is facial wrinkle formation. In recent years, an increasing number of patients have required more minimally invasive or nonsurgical procedures than surgical procedures for facial rejuvenation.^{1,2} Unwanted wrinkles can be treated by injection of dermal fillers, which fulfills the requirement of patients for minimal down time and which is a noninvasive treatment. Hyaluronic acid dermal fillers are the most frequently used class of injectable products for soft-tissue augmentation.³ Hyaluronic acid is a naturally occurring polysaccharide, consisting of linear chains of alternating D-glucuronic acid and N-acetyl-D-glucosamine

Disclosures: The Neuramis Deep Lidocaine and Restylane Perlane-L used in this clinical trial were generously supported by Medytox Inc., Korea. However, the authors have no financial interest to disclose in relation to the content of this article.

From the Department of Dermatology, St. Paul Hospital, College of Medicine, Catholic University of Korea; and the Department of Dermatology, College of Medicine, Chung-Ang University.
Received for publication May 11, 2015; accepted November 5, 2015.
This randomized, multicenter, double-blind, intraindividual, clinical trial was approved by the Ministry of Food and Drug Safety in the Republic of Korea. The registration number of MFTS is NCT02211348. This clinical trial was conducted in the Catholic Medical Center Clinical Research Coordination Center according to the guideline of the institutional review board of two centers, St. Paul's Hospital and Chung-Ang University Hospital in Seoul, Republic of Korea.
Copyright © 2016 by the American Society of Plastic Surgeons
DOI: 10.1097/01.prs.0000479965.14775.f0

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APS
Archives of Plastic Surgery

A Phase III, Randomized, Multi-Center, Double-Masked, Matched-Pairs, Active-Controlled Trial to Compare the Efficacy and Safety between Neuramis Deep and Restylane in the Correction of Nasolabial Folds

Changsik Pak¹, Jihoon Park¹, Jinmyung Hong¹, Jaehoon Jeong¹, Saik Bang², Chan Yeong Heo³

¹Department of Plastic and Reconstructive Surgery, Seoul National University Bundang Hospital, Seoul National University College of Medicine, Seongnam; ²Department of Plastic and Reconstructive Surgery, Samsung Medical Center, Sungkyunkwan University School of Medicine, Seoul, Korea

Background: We conducted this clinical study to compare the efficacy and safety between Neuramis Deep and Restylane in the correction of nasolabial folds.

Methods: In this phase III, randomized, multi-center, double-masked, matched-pairs, active-controlled trial (ClinicalTrials.gov Identifier: NCT01585220), we evaluated a total of 67 subjects (n = 67). All the subjects underwent Neuramis Deep treatment on one side and Restylane on the contralateral side of the bilateral nasolabial folds at a ratio of 1:1. To compare the efficacy of Neuramis Deep and Restylane, we evaluated the Wrinkle Severity Rating Scale scores and those of the Global Aesthetic Improvement Scale. In addition, we compared the safety of Neuramis Deep and Restylane based on adverse events, physical examination, and clinical laboratory tests.

Results: Neuramis Deep was not inferior in improving the nasolabial folds as compared with Restylane. In addition, there was no significant difference in the efficacy between Neuramis Deep and Restylane. There were no significant differences in safety parameters between Neuramis Deep and Restylane.

Conclusions: In conclusion, our results indicate that Neuramis Deep may be a safe, effective material for improving the nasolabial folds. However, further studies are warranted to compare the tolerability of Neuramis Deep and Restylane based on histopathologic findings.

Keywords: Hyaluronic acid / Nasolabial fold

Received: 3 Mar 2015 • Revised: 19 Jun 2015 • Accepted: 22 Jun 2015
ISSN: 2234-6163 • eISSN: 2234-6171 • http://dx.doi.org/10.1097/apc.0b013e3182e2721f • Arch Plast Surg 2015;42:721-728

INTRODUCTION

Hyaluronic acid (HA) is a polyanionic macropolysaccharide that is widely present in nature and extensively distributed in the connective tissues, including the skin and synovial fluid of the joints. It has high viscosity and hydrophilic properties, and is involved in many biological functions of the human body. However, it is promptly degraded by hyaluronidase, with a short half-life of

Correspondence: Chan Yeong Heo, Department of Plastic and Reconstructive Surgery, Seoul National University Bundang Hospital, Seoul National University College of Medicine, 82 Gumi-ro 173beon-gil, Bundang-gu, Seongnam 13620, Korea
Tel: +82-31-787-7222
Fax: +82-31-787-4055
E-mail: lonheo@bnuh.ac.kr

The current study was financially supported by a grant of the Korea Healthcare Technology R&D Project, Ministry of Health & Welfare, Republic of Korea (A100064).

No potential conflict of interest relevant to this article was reported.

Copyright © 2015 The Korean Society of Plastic and Reconstructive Surgery
This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (<http://creativecommons.org/licenses/by-nc/4.0/>) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited. | www.e-apj.org

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Received: 28 April 2019 | Revised: 25 September 2019 | Accepted: 29 October 2019
DOI: 10.1111/jcod.12018

ORIGINAL CONTRIBUTION

"A Multi-center, randomized, double blinded, comparative study of two hyaluronic acid fillers for temporary restoration of mid-face volume in Asians"

Joon Min Jung MD¹ | Woo Shun Lee MD² | Hak Tae Kim MD¹ | Young Seon Moon RN¹ | Seung Hwan Paik MD¹ | Woo Jin Lee MD, PhD¹ | Sung Eun Chang MD, PhD¹ | Mi Woo Lee MD, PhD¹ | Chong Hyun Won MD, PhD¹ | Beom Joon Kim MD, PhD^{2*}

¹Department of Dermatology, Asan Medical Center, University of Ulsan College of Medicine, Seoul, Korea
²Medytox Inc., Seoul, Korea
³Department of Dermatology, Chung-Ang University College of Medicine, Seoul, Korea

Correspondence: Beom Joon Kim, Department of Dermatology, Chung-Ang University College of Medicine, 84 Heukseok-ro, Dongjak-gu, Seoul 06974, Korea.
Email: beomjoonkim@ca.ac.kr

Chong Hyun Won, Department of Dermatology, Asan Medical Center, University of Ulsan College of Medicine, 88, Chongno-ro 43-gil, Songpa-gu, Seoul 05505, Korea.
Email: drwon@amc.seoul.kr

Funding information: Medytox Inc., Seoul, Korea

Abstract

Background: There have been few studies using hyaluronic acid fillers to restore mid-face volume deficit especially in Asians.

Aims: This study compared the efficacy and safety of two highly cohesive hyaluronic acid fillers, Neuramis® Volume Lidocaine and Juvederm® Voluma® with Lidocaine (VYC-20L), for temporary restoration of mid-face volume in Asians.

Patients/Methods: A total of 88 subjects with moderate to severe age-related mid-face volume deficit on the Mid-Face Volume Deficit Scale (MFVDS) received Neuramis® Volume Lidocaine on one side and VYC-20L on the other side of the face. Response was defined as a 1 point reduction in MFVDS at 24 weeks after the injection by photographic assessment. Secondary end points included photographic and live assessments of MFVDS and its changes at week 4, 12, and 24, global aesthetic improvements, as assessed by the investigators and the subjects at each visit, and the proportion of subjects who received touch-up treatment.

Results: The response rates of the sides treated with Neuramis® Volume Lidocaine and VYC-20L were identical (95.39%) at week 24, demonstrating that Neuramis® Volume Lidocaine was noninferior to VYC-20L. Photographic and live assessments of MFVDS and its changes, global aesthetic improvement, and the proportion of subjects who received touch-up treatment did not differ significantly during follow-up.

Conclusions: Neuramis® Volume Lidocaine was not inferior to VYC-20L in temporarily restoring mid-face volume at 24 weeks after treatment. Both of these highly cohesive hyaluronic acid fillers can be used effectively and safely for the correction of mid-face volume loss in Asians.

KEYWORDS
Asian, filler, hyaluronic acid, mid-face, Neuramis, volume

*These authors contributed equally as corresponding authors.
The study protocol was approved by the relevant institutional review boards and performed in compliance with the principles of the Declaration of Helsinki.
ClinicalTrials.gov Identifier: NCT02211348. URL: <https://clinicaltrials.gov/ct2/show/NCT02211348>

J Codemest. 2019;001-8. | wileyonlinelibrary.com/journal/jcod | © 2019 Wiley Periodicals, Inc. | 1

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